



Coccidioidomycosis – What’s the Best Test?

At a Glance – Key Points

1. Coccidioidomycosis is enzootic in dogs in the southwestern and western U.S.
2. Diagnosis is often based on clinical findings and a positive antibody test
3. *Coccidioides* antibodies can be detected via MVista® *Coccidioides* Canine IgG EIA and fungal immunodiffusion (FID)
4. The MVista® *Coccidioides* IgG Antibody EIA has the advantages more rapid turnaround time and being fully quantitative
5. Testing for antibodies with both the MVista® *Cocci* IgG EIA and FID antibody test is recommended for dogs suspected of having coccidioidomycosis

Case presentation

Presley, a 10-year-old spayed, female chihuahua from Arizona was presented for lethargy, fever (105F), dry cough, and dyspnea. The dyspnea was characterized by increased expiratory effort and an audible wheeze.



Figure 1: Transverse plane thoracic CT showing mainstem bronchial compression by enlarged tracheobronchial lymph nodes (calipers).

Thoracic radiographs showed probable tracheobronchial lymphadenopathy, which was confirmed on a thoracic CT. (Figure 1) Tracheoscopy showed extramural compression primarily of the main stem bronchi.

Coccidioidomycosis was suspected based on the clinical findings and that fact that Presley resided in an enzootic area. A *Coccidioides* antibody immunodiffusion (FID) test was submitted to a commercial laboratory and was negative. An additional FID test 4 weeks later was also negative. Due to the continued high suspicion for coccidioidomycosis, fluconazole (20 mg/kg/day) was prescribed. Over the following 12

months Presley’s clinical signs and lymphadenopathy resolved. (Figure 2) Follow-up testing of a stored serum sample from the time of diagnosis with a commercially available *Coccidioides* IgG antibody test (MVista® *Coccidioides* Canine IgG Antibody EIA) was positive (21.1 EU; >10 EU = positive).



Figure 2: Right lateral thoracic radiograph showing resolution of tracheobronchial lymphadenopathy after fluconazole treatment.

In Presley’s case, testing for antibodies with the MVista® *Cocci* IgG EIA at the time of presentation would have confirmed the diagnosis. Being quantitative it would have also provided a potential marker for treatment monitoring. Presley’s case highlights the recommendations discussed below, which include testing for *Coccidioides* antibodies with the MVista® *Cocci* IgG EIA and FID in dogs with suspected coccidioidomycosis.



CLINICAL DIAGNOSIS

Additional Information

Coccidioidomycosis is the most common invasive fungal infection (IFI) in dogs in the arid southwestern U.S., also being found in the western and southcentral U.S. Clinical signs are dependent upon the body systems affected, commonly limited to the lung and associated lymphatics causing cough, tachypnea, and dyspnea. Other common signs include fever, weight loss, anorexia, lameness and joint effusion, ulcerative or draining skin lesions, and seizures [1].

**MVista® *Coccidioides* IgG EIA:
low cost, quick
turnaround, quantitative**

Coccidioidomycosis is often diagnosed with appropriate clinical findings and a positive antibody test. Immunodiffusion is the traditional test for anti-*Coccidioides* antibodies (aka. agar gel immunodiffusion-AGID or fungal immunodiffusion-FID). The FID test detects precipitation when 2 *Coccidioides* antigens (test reagents) come into contact with antibodies (from the patient's serum) after migration through an agarose gel. The interpretation is somewhat subjective and can be delayed for up to 3 days after test set-up. Moreover, for positive results, titers require dilution of the sample and repeated testing requiring up to 3 additional days.

**MVista® *Coccidioides* IgG EIA and
FID testing should be done in dogs
suspected of having coccidioidomycosis**

More recently, MiraVista Diagnostics developed an immunoassay to detect anti-*Coccidioides* IgG antibodies, (MVista® Canine *Coccidioides* IgG Antibody EIA). This test improved upon the FID test with a lower cost, quicker result from time of test set-up (<12 hours), and being fully quantitative - all important considerations for diagnosis and treatment monitoring [2]. For example, the MVista® *Cocci* IgG EIA costs 32-45% less than the typical cost for the FID antibody test. Being quantitative, IgG antibody concentrations >10 EU are considered positive.

The *Coccidioides* antibody FID test gives a false negative result in 13% of dogs with biopsy confirmed coccidioidomycosis [3]. In order to diagnose those dogs, a recent study showed the combination of MVista® *Cocci* IgG EIA and antibody FID testing had a sensitivity of 98.6%, which was the highest combination in the study. These findings suggest that both tests should be performed in dogs suspected of having coccidioidomycosis. As with any antibody test, a positive result has to be interpreted considering the full clinical picture, as past exposure or asymptomatic infection can lead to the production of antibodies without necessarily requiring antifungal treatment.

REFERENCES:

1. Sykes J. Coccidioidomycosis. In: Sykes J, ed. *Canine and Feline Infectious Diseases*. St. Louis, MO: Elsevier; 2014:613-623.
2. Holbrook ED, Greene RT, Rubin SI, et al. *Novel canine anti-Coccidioides immunoglobulin G enzyme immunoassay aids in diagnosis of coccidioidomycosis in dogs*. *Medical mycology* 2019;57:800-806.
3. Gunstra A, Steurer JA, Seibert RL, et al. *Sensitivity of Serologic Testing for Dogs Diagnosed with Coccidioidomycosis on Histology: 52 Cases (2012-2013)*. *J Am Anim Hosp Assoc* 2019;55:238-242.

HEADQUARTERS

4705 Decatur Blvd. | Indianapolis, Indiana 46241, USA

888-841-8387