



## **Blastomycosis: When to Stop Antifungal Treatment.**

### **At a Glance – Crucial Information**

1. MVista® *Blasto* Antigen test can be used for diagnosis, treatment monitoring, and detection of disease relapse [1-4].
2. Treatment monitoring includes following clinical signs, physical exam findings, imaging study results, and antigen concentrations.
3. Antigen concentrations decrease with effective treatment and increase with disease relapse [4].
4. MVista® *Blasto* Antigen test should be performed at diagnosis, every 3 months during treatment, and at 6 months then every 12 months after treatment.

The MVista® *Blastomyces* Antigen test has revolutionized diagnosis and treatment monitoring of blastomycosis. It provides a non-invasive option with high diagnostic performance (sensitivity = 93%, specificity = 98%) [1-4]. With effective treatment, antigen concentrations decrease. With ineffective treatment or disease relapse, antigen concentrations remain static or increase over time. Many dogs and most cats have no detectable antigen at the time of disease remission. Less frequently, a very low concentration of antigen is present, even after remission. In these cases, it might be appropriate to stop treatment even when the antigen test is still positive (described below).

Urine is the ideal sample when using the MVista® *Blasto* Antigen test for treatment monitoring. It commonly contains a higher antigen concentration, as compared with serum, and thus is a better marker for remission and relapse. When urine antigen concentration is too high to be quantified (above the limit of quantification), serum can often be used instead until urine concentrations are quantifiable. In these cases, antigen is often quantifiable in urine at the 6 months recheck.

Disease relapse has been reported in approximately 20% of dogs treated for blastomycosis [5-8]. Intuitively, a treatment duration that is too short increases the chance of disease relapse. Using a multimodal monitoring approach helps identify the appropriate treatment duration (Table 1). In addition to treatment duration, appropriate drug and dosage are also vital. Itraconazole is the first-choice treatment for animals not hospitalized. Due to variability in absorption, itraconazole blood levels should be monitored to determine the ideal individualized dose. This can be achieved with the MVista® Itraconazole Bioassay, which accounts for parent drug and all active metabolites. It has the additional advantage of being significantly less expensive than chromatography/spectrometry testing available from pharmacology labs.

### **HEADQUARTERS**

4705 Decatur Blvd. | Indianapolis, Indiana 46241, USA

**888-841-8387**



## CLINICAL TREATMENT & MONITORING

**Table 1.** All of the following criteria should be met before stopping antifungal treatment in dogs and cats with blastomycosis.

Monitoring Tool	Criteria	Notes
<b>Treatment Duration</b>	Minimum of 6 months	<ul style="list-style-type: none"> <li>Required duration is often much longer.</li> </ul>
<b>History</b>	≥1-month past resolution of clinical signs	<ul style="list-style-type: none"> <li>Mild exercise intolerance might persist, most notable in working or performance animals.</li> <li>Persistent tracheobronchial lymphadenopathy can cause cough, requiring concurrent corticosteroid treatment.</li> </ul>
<b>Physical Examination</b>	≥1-month past resolution of physical exam abnormalities	<ul style="list-style-type: none"> <li>Differentiating active ocular disease from permanent inactive change is important.</li> </ul>
<b>Imaging Studies</b>	≥1-month past resolution of imaging abnormalities	<ul style="list-style-type: none"> <li>Pulmonary scarring can be permanent and can cause static focal unstructured interstitial lung disease.</li> <li>Radiographic bone lesions should improve but might never return to normal.</li> </ul>
<b>MVista® <i>Blasto</i> Antigen EIA (urine)</b>	≥1-month past no detectable antigen OR Antigen ≤0.4 ng/ml on 2 consecutive rechecks at least 3 months apart	<ul style="list-style-type: none"> <li>Most dogs, and essentially all cats, have no detectable antigen at the time of remission.</li> <li>Submit antigen test at diagnosis, every 3 months during treatment, and at 6 months then every 12 months after treatment.</li> </ul>

### REFERENCES:

- Motschenbacher LO, Furrow E, Rendahl AK, et al. Retrospective analysis of the effects of *Blastomyces* antigen concentration in urine and radiographic findings on survival in dogs with blastomycosis. *J Vet Intern Med* 2021.
- Mourning AC, Patterson EE, Kirsch EJ, et al. Evaluation of an enzyme immunoassay for antibodies to a recombinant *Blastomyces* adhesin-1 repeat antigen as an aid in the diagnosis of blastomycosis in dogs. *J Am Vet Med Assoc* 2015;247:1133-1138.
- Spector D, Legendre AM, Wheat J, et al. Antigen and antibody testing for the diagnosis of blastomycosis in dogs. *J Vet Intern Med* 2008;22:839-843.
- Foy DS, Trepanier LA, Kirsch EJ, et al. Serum and urine *Blastomyces* antigen concentrations as markers of clinical remission in dogs treated for systemic blastomycosis. *J Vet Intern Med* 2014;28:305-310.
- Legendre AM, Rohrbach BW, Toal RL, et al. Treatment of blastomycosis with itraconazole in 112 dogs. *J Vet Intern Med* 1996;10:365-371.
- Arceneaux KA, Taboada J, Hosgood G. Blastomycosis in dogs: 115 cases (1980-1995). *J Am Vet Med Assoc* 1998;213:658-664.
- Mazepa AS, Trepanier LA, Foy DS. Retrospective comparison of the efficacy of fluconazole or itraconazole for the treatment of systemic blastomycosis in dogs. *J Vet Intern Med* 2011;25:440-445.
- Krawiec DR, McKiernan BC, Twardock AR, et al. Use of an amphotericin B lipid complex for treatment of blastomycosis in dogs. *J Am Vet Med Assoc*

### HEADQUARTERS

4705 Decatur Blvd. | Indianapolis, Indiana 46241, USA

888-841-8387